

# Methodological and analytical challenges for PFAS monitoring in environmental media

Update from the frontline

Stefan Voorspoels, Griet Jacobs



## Agenda


- PFAS: more than 1 compound (class)
- PFAS: more than 1 analytical method
- Challenges
- Possible solutions
- Caveats and restrictions
- Answers or questions ?

*“you can only find what you search for”*

# PFAS: More than 1 compound (class)


## WHAT?

- Molecules with varying (i) carbon chain length, (ii) functional groups and (iii) number of fluor atoms
- PFAA + Polymers + Precursors + **polyfluorinated cpds**
- Water, soil, air, **soil improver, waste, consumer products, industrial intermediates, animals, food, humans, ...**
- Particle-bound, semi-volatile and volatile



<https://doi.org/10.3390/toxics10020044>

## Per- and Polyfluoroalkyl Substances (PFAS; $C_nF_{2n+1}-R$ )




<https://doi.org/10.1029/2021RG000765>



vito.be


# PFAS: More than 1 compound ... more than 1 classification

OECD 2020 classification



Pieter Bruegel de Oude - [https://en.wikipedia.org/wiki/The\\_Tower\\_of\\_Babel\\_\(Bruegel\)](https://en.wikipedia.org/wiki/The_Tower_of_Babel_(Bruegel))

## Per- and polyfluoroalkyl substances (PFASs)




# PFAS: More than 1 compound ... more than 1 classification




Pieter Bruegel de Oude -  
[https://en.wikipedia.org/wiki/The\\_Tower\\_of\\_Babel\\_\(Bruegel\)](https://en.wikipedia.org/wiki/The_Tower_of_Babel_(Bruegel))

- Unique ID
- IUPAC
- PFOS = 1,1,2,2,3,3,4,4,5,5,6,6,7,7,8,8,8-heptadecafluorooctane-1-sulfonic acid
- CAS #
- e.g. ADONA, NaDONA, DONA
- What with:
  - Polymers
  - Break-down products
  - Derivatives
  - PolyF-cpds
  - ...
- Quantification & Reference compounds ???




<https://www.facebook.com/groups/encyclopedia.nonsensical/>

# PFAS: More than 1 compound ... more than 1 classification




PFBA-  
precursor

Chemical Formula:  $C_{20}H_{21}F_{14}O_{11}S^-$   
Exact Mass: 735,0587




Chemical Formula:  $C_{15}H_{17}F_{13}NO_6S_2^-$   
Exact Mass: 618,0295

Chemical Formula:  $C_{27}H_{28}F_{21}O_{12}S^-$   
Exact Mass: 975,0972



Chemical Formula:  $C_{11}H_{10}F_{13}O_4S^-$   
Exact Mass: 485,0098  
Acetate adduct



Chemical Formula:  $C_{27}H_{28}F_{21}O_{12}S^-$   
Exact Mass: 975,0972

# PFAS: More than 1 method

Intended for  
NORDIC WORKING GROUP FOR CHEMICALS, ENVIRONMENT, AND HEALTH (NKE)

Document type  
Final Report

Date  
December 2023

## ANALYSIS OF NEEDS FOR ENFORCEMENT OF PFAS IN ARTICLES AND CHEMICAL PRODUCTS



© chakniti - stock.adobe.com

### Section 2.1

#### Total fluorine content

- CIC
- LIBS
- WDXRF
- PIGE
- INAA
- HR-CS-GF-MAS
- XPS

All fluorinated compounds  
No individual identification/quantification

### Section 2.2

#### Non-targeted and suspect screening analysis via HRMS

- LC- and GC-MS(MS)
- MALDI
- IM-MS
- SPE-coupling
- DART
- ICP-MS
- direct infusion ESI
- FT-ICR

Several PFAS simultaneously  
High sensitivity and specificity

### Section 2.3

#### Targeted analysis

- LC- and GC-MS
- SFC-MS
- Sensors
- pyr-GC-MS
- Sensors
- TOPA

Individual PFAS  
High sensitivity and specificity

### Resource quantifiers

- Time investment
- Costs
- Expertise

### Section 2.4

#### Structural analysis & additional methods

- <sup>19</sup>F-NMR
- XANES
- FTIR
- MIP-OES
- SERS

Individual PFAS  
High specificity

→ quantification possible  
... if reference compounds available


# PFAS: More than 1 method

## Quantification of PFAS – Target analysis

- Multiple methods needed: legacy cpds, ultra-short chain cpds, volatile PFAS, ... !
- It will **NEVER** be possible to quantify all PFAS we encounter
- Applicability of target method is different according to corresponding matrix and PFAS class
  - WAC/IV/A/025 – Determination of PFAS in water with LC-MS/MS - #43
  - WAC/IV/A/026 – Determination of ultra-short chain PFAS in water with LC-MS/MS - #7
  - CMA/3/D – Determination of PFAS in soil and sediment - #43
  - CMA/3/O – Determination of PFAS in soil improvers - #42
- Expected updates
  - Suspect and non target analysis (S/NTA) → selection of new relevant target compounds
  - Depends on the availability of reference standards and internal standards



## PFAS: More than 1 method



Romolo Tavani | Dreamstime.com

# Challenges and solutions

## Total PFAS methods


### The GOOD

- Measuring total F, but for EOF and AOF the inorganic fluorine should be removed (extraction/adsorption)
- Can give a quick idea of the total fluorine content
- ISO standard is under development (AOF)

### The BAD

- Not all techniques distinguish between inorganic and organic F (false positives)
- The limit of F detection mg/L range
- Molecular technique vs SI-data reporting
- Not (yet) wide-spread in commercial labs.
- CIC and HR-CS-GF-MAS can be used but these are destructive methods where often sample preparation is needed
- Ultrashort-chain PFAS remain a blind spot even for these “total PFAS” methods

## TOF-techniques



[Gehrenkemper et al., 2021](#)

[DOI: 10.1007/s00216-020-03010-y](https://doi.org/10.1007/s00216-020-03010-y)

# Challenges and solutions

## Suspect and non target analysis (S/NTA)


### The GOOD

- HRMS allows identification of new and suspected PFAS
- HRMS can measure both multiple known and unknown analytes in a single analysis
- LOD factor 10 or more lower than target analysis
- Developments for high-throughput ongoing
- Results will direct future developments of target methods

### The BAD

- High cost of technology acquisition, operation and maintenance
- High training and experience level of scientists
- Large datasets
- Slower process than target analysis
- Not (yet) assimilated in commercial labs

## NTA-techniques




[Charbonnet et al., 2022](#)

[DOI: 10.1021/acs.estlett.2c00206](https://doi.org/10.1021/acs.estlett.2c00206)

# Challenges and solutions

- Complementary techniques are needed
- Orthogonal approaches
- Complexity tailored to the need
- Target methods will not become obsolete but will require maintenance and updates
- ... the inherent limitations of every analysis applies ...



# Challenges and solutions

Update from the frontline

- Expansion of Targeted methods LC-MS/MS
  - More compound classes
- Development, validation and implementation of TOP assay
- Development/improvement of GC-MS method for the analysis of neutral and volatile PFAS
- Development/improvement of AOF and EOF (total PFAS methods)
- Further development and application of the NTA PFAS platform
- Development of a methodology for the quantification of PFAS in flue gasses
- ...





**Uncertainty**  
Just Ahead

***“Uncertainty is  
unpleasant, but  
certainty is absurd”***

Voltaire (1694-1778)




<https://en.wikipedia.org/wiki/Voltaire>





# Caveats and restrictions


- Uncertainty of:
  - Compound ID
  - Compound quantity



Bijlage bij accreditatiecertificaat  
Annexe au certificat d'accréditation  
Annex to the accreditation certificate  
Beilage zur Akkreditierungszertifikat

**045-TEST**  
EN ISO/IEC 17025:2017

## fit-for-purpose measurement uncertainty



**Meta-data**

**5 %**  
**20 %**  
**50 %**






**Analytical result:**

$$x \pm U$$

vito.be

# Caveats and restrictions – uncertainty

sampling – preparation – measurement – data interpretation




## Caveats and restrictions – uncertainty



# Caveats and restrictions – uncertainty

## The curious case of PFOS-isomers (aka branched vs linear)



NOTE: R= CO<sub>2</sub><sup>-</sup> and CF<sub>2</sub>SO<sub>3</sub><sup>-</sup>  
R'= SO<sub>3</sub><sup>-</sup> only

# Caveats and restrictions – uncertainty

## The curious case of PFOS-isomers (aka branched vs linear)




## Caveats and restrictions

- Knowing what to count → > 6000 PFAS
- Specificity matters
- **Isomers – linear vs branched**
- Qualitative analysis → identity + present or not (LOD)
  - How to identify?
- Quantitative analysis → identity + mass fraction (amount)
  - How to identify?
  - How to quantify? → calibration – **reference (SI-traceable)** – reproducible data (comparison)

## Caveats and restrictions

- Reference compounds of known purity needed
  - huge response factor differences, even for closely correlated structures
- For true quantification: matching mass-labelled IS standards essential
- Every analytical measurement is only an estimation
  - Measurement uncertainties range from 15 – 50 % ( $k = 2$ )
  - But: Fit-for-purpose data
- Setting realistic expectations (read: quality)



***“Analytical chemistry is an exact science,  
as long as it is exact enough”***

[stefan.voorspoels@vito.be](mailto:stefan.voorspoels@vito.be)

*Disclaimer: The information in this communication is the result of professional research and is based on advancing insights and empirical knowledge of high-quality scientific and technical research by VITO. Since this information is based on the state of VITO's scientific, technical and policy knowledge at the time of publication, the completeness and accuracy of this information cannot be guaranteed for future developments. The research into PFAS is complex; We guarantee the accuracy of our figures, but we cannot accept responsibility for conclusions that third parties base on these figures.*



**vito.be**